


# Manually installing USB drivers


## (MA60I/MA80I connected to Charder Insight) (Windows)

- Scenario:** Charder Insight is unable to connect with MA60I/MA80I Body Composition Analyzer, even though UA-010 USB Transfer Key Cable is already plugged in. USB driver needs to be installed manually on PC.
- Applicable operating systems:**  
Windows 10, Windows 11
- Checking if drivers need to be installed:**
  - Plug UA-010 Cable into Windows PC
  - Press [Windows Key] + [R] to open the “Run Command”. Type “devmgmt.msc” and press [Enter].
  - Under “Other devices”, confirm if an exclamation mark “!” is seen next to “USB Serial Port”. If the exclamation mark appears, then USB drivers need to be installed.




#### 4 Installing USB driver:

- After downloading CDM212364\_Setup.zip file from [Charder Medical website](#), extract “CDM212364\_Setup.exe” and run installation program.



- Follow the installation instructions:.



Device Driver Installation Wizard



## Welcome to the Device Driver Installation Wizard!



This wizard helps you install the software drivers that some computers devices need in order to work.

To continue, click Next.

< Back **Next >** Cancel

Device Driver Installation Wizard

### License Agreement



To continue, accept the following license agreement. To read the entire agreement, use the scroll bar or press the Page Down key.

**IMPORTANT NOTICE: PLEASE READ CAREFULLY BEFORE INSTALLING THE RELEVANT SOFTWARE:**  
This licence agreement (Licence) is a legal agreement between you (Licensee or you) and Future Technology Devices International Limited of 2 Seaward Place, Centurion Business Park, Glasgow G41 1HH, Scotland (UK Company Number SC136640) (Licensor or we) for use of driver software provided by the Licensor(Software).

BY INSTALLING OR USING THIS SOFTWARE YOU AGREE TO THE

**accept this agreement!**  I don't accept this agreement

Save As Print

< Back **Next >** Cancel



- After driver is installed, exclamation mark “!” will disappear. Confirm if MA60I/MA80I can connect to Charder Insight now.

