



CUPID 5 Digital Infant Scale

Features

- Easy-to-read large 1-inch LCD display
- Tare feature that accounts for extra weight attributed to diaper, clothing or blanket
- It includes power saving feature ease your pain of replacing batteries
- Sleek contour design makes baby feel safe and comfortable
- Blue backlight display, ensure clearly measurement even in dark
- High accuracy weighing sensor ensures every measurement

Specifications

Capacity	20 kg / 44 lb
Graduation	0.01 kg / 0.02 lb
Measuring unit	kg and lb
Weight	1.6 kg
Dimensions	55x33x4 cm
Power supply	4 x AA batteries
Display type	LCD with backlit



Charder Electronic Co., Ltd.

No. 103, Guozhong Rd., Dali Dist., Taichung City 41262, Taiwan

Tel 886-4-2406-3766

Fax 886-4-2406-5612

info_cec@charder.com.tw

www.chardermedical.com