



Carrying case



Wireless Grip Strength Dynamometer

The new All-Age variant of the MG4800 Wireless Grip Strength Dynamometer is specially designed to be used by both children and adults for precision grip strength measurement. When used by adults that require a larger grip span, a grip extension can be added in seconds for comfortable usage!

Three modes are available, measuring maximum grip strength, average grip strength, or number of repetitions that meet the selected minimum grip force.

Results can be transmitted wirelessly to EMR system for more convenient data management and tracking.

Featuring multi-color charging status indicator for rechargeable battery, a full charge can be completed in 2.5 hours.

Model	MG4800 (All-Age)
Capacity	80 kg
Grip span	64 mm
Graduation	0.1 kg
Overall Dimensions	110 (L) x 212 (D) x 55 (W) mm
Device Weight	0.33 kg
Power Supply	CR123A Li-ion rechargeable battery
Data Transfer	Wireless & USB
Accessories	Carry case, USB Type-C cable, USB charger, Grip extension
Functions	<ul style="list-style-type: none"> • Average grip strength • Maximum grip strength • Count Mode