




Versatile Needs

MHS2500I

Lift Scale



Designed to meet different lift/mobility systems



Features

- Large size LCD display
- Hold function to lock weight of moving patient
- Joint protection design
- Can be used with most manufacturers lifts
- Compact size and reliable performance

Specifications

Model	Display	Capacity	Graduation	Dimensions	Functions	Data Transmission
MHS2500I	Embedded 1.3 inch LCD screen	300kg 400kg	100g 200g	52 x 122 x 180 mm	Power on/off, zero, UNIT, Hold, BMI, tare	N/A