

# Digital Veterinary Scale MS2200



## Features

- 300 kilograms (660 pounds) capacity
- Measures in kilograms or pounds
- Designed with high-precision load cells
- Large LCD display with backlit
- Easy-to-clean stainless steel platform and removable non-slip mat
- Hold function locks weight measurements of a moving animal
- Tare function compensates for weight of owner or pad
- Scale powered by 4 AA batteries(not included) or AC adaptor(included)
- Carry handle and caster wheels for portability



Removable, Washable Mat

## Specifications

|            |                      |
|------------|----------------------|
| Capacity   | 300 kilograms        |
| Graduation | 0.1 kilograms        |
| Dimensions | 90 x 55 x 7 cm LxWxH |
| Accuracy   | +/- 0.2 kilograms    |
| Weight     | 21 kilograms         |