

Digital Veterinary Scale MS2210R



Features

- 300 and 150 kilograms capacity
- Measures in kilograms or pounds
- Designed with high-precision load cells
- Remote wireless control panel with large LCD display
- Easy-to-clean stainless steel platform and removable non-slip mat
- Hold function locks weight measurements of a moving animal
- Tare function compensates for weight of owner or pad
- Scale powered by 4 AA batteries(not included) or AC adaptor(included)
- Carry handle and caster wheels for portability



Removable, Washable Mat

Specifications

Capacity	300 or 150 kilograms
Graduation	0.1 or 0.05 kilograms
Dimensions	90 x 55 x 7 cm LxWxH
Accuracy	+/- 0.2 or +/- 0.1 kilograms
Weight	21 kilograms