

MS3400-1 + HM201M

Digital Column Scale with Mechanical Height Measurement

Minimalistic and compact, the MS3400-1 features a waist-level column, large LCD screen, and steel platform with anti-slip coating to make it an economic medical scale for busy clinics.

Optional height stadiometer can be attached to make easy use of the scale's integrated BMI function.

Model with optional castor wheels also available!

Functions

- BMI
- Unit (kg / lb)
- Hold
- Rechargeable Battery
- Tare
- Castor for mobility (optional)
- Zero

Specifications

Capacity	300 kg
Graduation	0.1 kg
Overall Dimensions	460(D) x 310(W) x 1070(H) mm (Standard Model)
	504(D) x 310(W) x 1260(H) mm (Castor Model w/Stadiometer Fully Retracted)
Device Weight	5.3 kg (Standard Model) 6.6 kg (w/Stadiometer & Castor)
Power Supply	Rechargeable Battery / AC Adapter
Height Measurement	HM201M: Manual Mechanical Stadiometer HM200D: Manual Digital Stadiometer (LCD) HM201D: Manual Digital Stadiometer
Data Transmission	RS232

HM201M Mechanical Height Stadiometer

Capacity	60-210 cm
Graduation	0.1 cm



Foldable head stopper

Measures weight and BMI

Castor Wheel (optional)

Image contains optional
Mechanical Stadiometer (HM201M)