

MS3450 + HM201M

Digital Column Scale with Mechanical Height Measurement

Minimalistic and compact, the MS3450 features a waist-level column, large LCD screen, and steel platform with anti-slip mat to make it an economic medical scale for busy clinics.

Optional height stadiometer can be attached to make easy use of the scale's integrated BMI and BSA functions.

Model with optional castor wheels also available!

Functions


- BMI
- BSA
- Print
- Hold
- Tare
- Unit (kg / lb)
- Send
- Castor wheel (optional)

Specifications

Capacity	300 kg
Graduation	0.1 kg
Overall Dimensions	460(D) x 310(W) x 1070(H) mm (Standard Model) 504(D) x 310(W) x 1260(H) mm (Castor Model w/Stadiometer Fully Retracted)
Device Weight	5.3 kg (Standard Model) 6.6 kg (w/Stadiometer & Castor)
Power Supply	Rechargeable Battery / AC Adapter
Height Measurement	HM201M: Manual Mechanical Stadiometer HM200D: Manual Digital Stadiometer (LCD) HM201D: Manual Digital Stadiometer
Data Transmission	USB

HM201M Mechanical Height Stadiometer

Capacity	60-210 cm
Graduation	0.1 cm



Foldable head stopper



Measures weight, BMI, BSA



Castor wheel (optional)

Image contains optional Mechanical Stadiometer (HM201M)