




MS3500

Baby Scale with Removable Tray


The MS3500 infant digital scale features sloped edges to safely and comfortably hold infants during weighing. Highly accurate, the scale is able to track even slight changes in infants' weight. With a removable baby tray and 20kg capacity, it can be used to weigh toddlers as well.

Receive height results using the optional stadiometer, and send results wirelessly to compatible EMR systems!

View height results on LCD screen (and transfer results)
when compatible infant stadiometers (HM80DT, HM100D) are attached



Transmit results wirelessly to Charder Connect



Functions

- Hold
- Tare
- Removable Tray
- Wireless Data Transfer(optional)

Specifications

Capacity	20 kg
Graduation	5g < 10 kg > 10g
Overall Dimensions	560(W) x 342(D) x 105(H) mm
Platform Dimensions	320(W) x 310(D) x 35(H) mm
Tray Dimensions	560(W) x 290(D) x 65(H) mm
Device Weight	2.6 kg
Accessories	Infant Height Stadiometer AR-00001 Carrying Bag
Data Transfer	Wireless (optional)
Software	Charder Connect

Option

Height Measurement	HM80D, HM80M ,
	HM80DT, HM100D (Bluetooth)



Measure toddlers with measurement tray removed



AR-2491
Optional carry bag for baby scale.