



MS4200

Infant Scale with Removable Tray

The MS4200 Digital Infant Scale features sloped edges to safely and comfortably hold infants during weighing. Highly accurate, the scale is able to track even slight changes in infants' weight. An optional infant height Stadiometer can be attached for height measurement, and weight data can be transferred wirelessly to printer. The measurement tray is easily removed, making it possible to weigh toddlers as well!

Functions

- Hold
- Tare
- Removable Tray
- Wireless Data Transfer(optional)

Specifications

| | |
|---------------------|--|
| Capacity | 20 kg |
| Graduation | 5g < 10 kg > 10g |
| Overall Dimensions | 560(W) x 325(D) x 145(H) mm |
| Platform Dimensions | 325(W) x 310(D) x 50(H) mm |
| Tray Dimensions | 560(W) x 290(D) x 65(H) mm |
| Device Weight | 4.3 kg |
| Accessories | Infant Height Stadiometer AR-249 I Carrying Bag |
| Data Transfer | Wireless (optional) |
| Software | Charder Connect |



AR-249 I
Optional carry bag for baby scale.



Compatible with HM80-series
Infant Height Stadiometers

Accessories

| | |
|--------------------|--------------|
| Height Measurement | HM80D, HM80M |
|--------------------|--------------|