



# MS5460

## Compact Chair Scale



Swivelling armrests



Foldable footrests



Compact size for delivery



Tool-free assembly design

### Features

- Four castor wheels with brakes
- Folding frames collapse for easy transportation
- One-person assembly, tool-free clip fastenings
- Able to output data to PC



#### Key Functions

- Hold
- Tare / Zero
- Body Surface Area
- Body Mass Index

### Specifications

Model	Display	Capacity	Graduation	Dimensions	Functions	Data Transmission
MS5460	DP4400 1.4 inch LCD screen	250 kg	100 g	590(H) x 1020(D) x 980(W) mm	BMI, BSA, Hold, Weight Measurement, Tare, Zero	USB