



MS5480

Wireless Compact Chair Scale

The MS5480 compact chair scale is specially designed for patients that have difficulty standing, allowing for the simple and precise weighing of patients. It is intended to be pushed directly into patient rooms, so patients do not need to walk to weighing stations. After completing a round of measurements (up to 100 stored), results can be transmitted wirelessly to PC, making the MS5480 an excellent option for integration with EMR systems.

Featuring collapsible folding frames for easy transportation, it can be assembled easily by one person without tools.



Swiveling armrests



Foldable footrests



Compact size for delivery

Functions

- BMI
- Pre-Tare
- Rechargeable battery for mobile use
- Hold
- Color Touchscreen
- Foldable frame
- Tare
- Castor Wheel brakes



Specifications

Capacity	250 kg
Graduation	0.1 kg
Overall Dimensions	990(D) x 590(W) x 960(H) mm
Seat Dimensions	Seat Height: 536 mm Seat Width: 460 mm
Device Weight	17.6 kg
Accessories	Barcode Scanner
Data Transfer	USB, Wireless
Software	Smart Data Manager