

# MS7301 Bariatric Floor Scale



## Features

- 550 lbs (250 kg) Capacity
- Automatically turns on and weighs when stepped on
- Automatic off
- Reliable, professional-grade measuring
- Wide platform is stable and accommodates bariatric patients
- Low profile allows for easy access
- Self adjusting feet provide prevent tipping
- Non-slip texture



MS7301 Top View

## Specifications

Dimensions	16" x 14" x 1.6" (415 x 360 x 41 mm) WxDxH
Graduation	0.2 lbs (100 g)
Power	(4) AA Batteries (Included)
Capacity	550 lbs (250 kg)