

## VS2400 Electronic Veterinary Scale

The Charder VS2400 is a professional and user-friendly veterinary scale which offers a simple way on veterinary care. This scale features easy-to-install that provides users with convenience and effortless when using. It also simplifies the complexity of services without any calibration after replacing the display. With a sided handle and an extremely light platform, it greatly increases the mobility that means users can weigh their beloved pets anywhere.





## Features :

- Large and easy-to-clean stainless steel platform
- Clear and large LCD display for easy reading
- Light and thin design for space saving
- Easy to operate 4 keys panel design
- Hold function locks weight reading when your pet is restless
- Tare function allows pets to be weighed with its owner with a mat, a box, or other objects



Ergonomic handle allows to carry easily



Plug and weigh function is design for install and replace the display easily.



Extremely thin platform that is less than 25 mm

## Specifications :

Model	VS2400
Capacity	180 kg
Graduation	100 g
Measuring unit	kg or lb
Weight	6.2 kg
Dimensions	Platform: 900 x 550 x 25 mm
Power supply	Adapter & 4 x AA Batteries
Key Functions	ON/OFF, HOLD, TARE, UNIT