

# CUPID 5 Digital Infant Scale






Max. 20kg



Division  
10 g / 0.02 lb



Ultra-slim



Tare  
function



Large LCD



Backlight

## CUPID 5 Digital Infant Scale

### FEATURES

- Easy-to-read large 1-inch LCD screen
- Tare feature makes it easy to deduct extra weight (ex: diaper, blanket)
- Power-saving auto-off feature
- Sleek curved contour design, making measurement safer
- Blue backlit display for improved legibility even in dim environments
- High accuracy weighing sensor ensures measurement accuracy

### SPECIFICATIONS

Capacity	20 kg / 44 lb
Graduation	0.01 kg / 0.02 lb
Measuring unit	kg and lb
Weight	1.6 kg
Dimensions	55 x 33 x 4 cm
Power supply	4 x AA batteries
Display type	Backlit LCD



### Charder Electronic Co.,Ltd.

NO.103, Guozhong Rd., Dali Dist., Taichung City 41262, Taiwan

Tel 886-4-2406-3766

Fax 886-4-2406-5612

info\_cec@charder.com.tw

www.chardermedical.com